Instructions for Use AT51113

Arthritisrollator AT51113

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Instrucciones de uso andador AT51113

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The AT51113 rollator is made of high-quality, durable aluminium. The folding frame makes transport easier. The rollator is equipped with front swivel castors and rear fixed castors with a brake and a parking clip. The product has an additional seat and backrest, which enables the user to rest in a sitting position. The rollator is also equipped with a convenient bag for essential items or small groceries.

Warnings

Check the condition of the device before performing any activities. If the product functions incorrectly or looks defective, please contact the dealer from whom it was purchased immediately. Antar assumes no responsibility for improper use of the device.

Product composition:

frame	plastic
aluminium	armrest padding
castors	polyurethane
PVC	seat
handles	nylon and polyester

Preparation for use

The type of aid device and its appropriate configuration must be selected by the qualified medical and technical personnel based on disease and method of using the device by the user.

1. Unfold the frame by lowering the transverse tubes of the seat.
2. Adjust the height of the handles and armrests.
3. Attach the bag.
4. Attach the centre belt (Fig. 1), make sure the belt tension is correct to prevent the rollator from folding accidentally.

Using the brakes


Fig. 1


Fig.2

1. To release or stop the rollator, pull on the brake lever. Release the brake lever to move again (Fig. 1).
2. To activate the parking brake, push the brake lever down until it locks (Fig. 2).

Maximum user's weight	136 kg
total width	61.5 cm
space between handles	50 cm
total depth	70 cm
seat width	46 cm
seat depth	25 cm
seat height	54 cm
handle height adjustment range	108.5 - 121 cm
castors	8"
weight	8.4 kg

THIS SIGN INDICATES THE MAXIMUM WEIGHT OF THE USER

Storage and disposal

Store the rollator away from sources of light, heat and humidity. Avoid exposing the rollator to direct sunlight. Dispose of the rollator in accordance with applicable environmental protection and recycling regulations.

General precautions during use

Before use, always make sure that the handles are at the correct height and are properly secured. To avoid the risk of falling, do not use the rollator on uneven or inclined surfaces. Check the condition of the castors and brakes. Finally, check the overall stability of the frame.

Maintenance and cleaning

After each use, remove accumulated dirt using mild detergents, then leave the rollator to dry. Do not use sharp rubbing sponges. Periodically check the operation of brakes, the condition of fastening elements, gaskets and welds.

WARNING: When using and operating the rollator and when folding and adjusting the mechanisms, there may be a danger of trapping and/or squeezing the user's / accompanying person's body parts in the holes/slots between the elements.

These operations should be performed with particular care.

Once adjustment is complete, stabilize the position by tightening nuts/screws carefully.

WARNING:

In the event of pain, allergic reactions or other alarming, unclear symptoms related to the use of the medical device, please consult a healthcare professional.

WARNING:

In the event of a device-related "serious incident", which directly or indirectly led to, could have led to, or is likely to lead to any of the following events:

- death of a patient, user or another person, or
- temporary or permanent deterioration of the health of a patient, user or another person, or
- a serious risk to public health, such a "serious incident" must be notified to the manufacturer and to the competent authority of the Member State, in which the user or patient resides. For Poland, the competent authority is the Office for Registration of Medicinal Products, Medical Devices and Biocidal Products.

WARNING: The manufacturer is not responsible for damage caused by negligent maintenance, inadequate servicing or failure to follow the guidelines contained in these Instructions for Use.

WARNING: The product must not be used outside its intended purpose!

WARNING: If the product is used incorrectly, there is a risk of tipping over. Please follow the guidelines for getting on/off/moving around. Once adjustment is complete, stabilize the position by tightening nuts/screws carefully.

WARNING: When using and operating the product and when folding and adjusting the mechanisms, there may be a danger of trapping and/or squeezing the user's / accompanying person's body parts in the holes/slots between the elements. These operations should be performed with particular care. Once adjustment is complete, stabilize the position by tightening nuts/screws carefully.

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Model:

Serial number:

Point of sale stamp
and a legible signature of the Seller:

1. ANTAR Sp. J., 03-068 Warsaw, ul. Zawiańska 43, hereby grants a 12-month warranty on the rollator, starting from the date of purchase by the Buyer. Products financed by the National Health Fund are covered with the extended warranty period, which covers half of the useful life specified in the Regulation on Guaranteed Services, provided that the useful life is specified in the Regulation in months or years.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, send it immediately to ANTAR (as specified in point 1) or contact the point of sale.
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum permissible weight of the user as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts subject to normal wear and tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED. IN ORDER TO EXTEND WARRANTY PERIOD FOR PRODUCTS FINANCED BY THE NATIONAL HEALTH FUND, ATTACH A COPY OF THE ORDER FOR THE SUPPLY OF MEDICAL DEVICES BEING ORTHOPAEDIC ITEMS AND AIDS.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

ANTAR Sp. J.
03-068 Warszawa, ul. Zawiańska 43
Polska
Tel. 22 518 36 00, Fax 22 518 36 30
www.antar.net e-mail: antar@antar.net

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Date of issue of the manual: April 1, 2022
v1-1 04/01/2022